

# Wasp

## BARCODE TECHNOLOGIES

Portal > Knowledgebase > Software Updates & Patches > WaspLabeler db Tracker > WaspLabeler db Tracker Updates & Patches

---

### WaspLabeler db Tracker Updates & Patches

Beau Dixon - 2021-02-12 - in WaspLabeler db Tracker

#### **NDC Database Update 1.0.0.18**

The FDA updated the NDC database on September 30th 2013. This download will update your WaspLabeler db Tracker NDC data to match the data from the FDA.

**Note:** This requires **WaspLabeler db Tracker 1.2.0.2**. If you have an older version, please install dbTrackerPatch.exe Version 1.2.0.2

File Name: **NDDFInstallV2.exe**

Version: **1.0.0.18**

Date Published: **10/01/2013**

Language: **English**

Download Size: **13.5 MB**

[Download Now](#)

### WaspLabeler db Tracker 1.0.0.16

The FDA released a new NDC database on June 6th 2012. This download will update your WaspLabeler db Tracker NDC data to match the data from the FDA.

**Note:** This requires **WaspLabeler db Tracker 1.2.0.2**. If you have an older version, please install **dbTrackerPatch.exe Version 1.2.0.2**

File Name: **NDDFInstallV2.exe**

Version: **1.0.0.16**

Date Published: **08/09/2012**

Language: **English**

Download Size: **6 MB**

[Download Now](#)

### WaspLabeler db Tracker 1.2.0.2

This updates WaspLabeler db Tracker to be compatible with the latest formats of the FDA release of the NDC database.

File Name: **dbTrackerPatch.exe**

Version: **1.2.0.2**

Date Published: **01/13/2012**

Language: **English**

Download Size: **6 MB**

[Download Now](#)

### WaspLabeler db Tracker 1.2.0.1

1.2.0.1 - June 2011

1. Fixed the "Repackaged Medications Report" to not display duplicate entries.
2. Fixed the "Approval Report" to show the approver name instead of showing the person who printed it.
3. Fixed the 64-bit compatibility. Now the WaspLabeler db Tracker is fully 64 bit compatible.
4. Updated the SQL Express installation to SQL Express 2005 SP4.
5. Included the latest version of Labeler 6.0.10 in the installer.

---

## **WaspLabeler db Tracker 1.2.0**

1.2.0 - October 21, 2008

1. Printing to label printers will now push the label size to the printer. This means you will not have to edit the printer driver settings every time you change media size in the printer. If the printer driver media size does not match the actual media size the print preview will be incorrect but the label should still print correctly.
2. The update link on the main page will now take you to the proper update page and download link when a new version is available.
3. In the documents folder on the CD is a file that describes each report and gives some tips for using them.
4. The approval and print history lists now show the user who logged in as First Name Last Name (used to be Last, First).
5. Template 7062A now has all the default settings correct.
6. User logins with restricted access will now be able to create and duplicate drugs. db Tracker must be run once with a user with full admin access.
7. Selecting Reprint from the print history or Approval list will now ask if you want to reject the print that is being reprinted and automatically re
8. ject it if you choose. If you do not have approvals turned on you will not be prompted for this.
9. The manufacturer search list no longer has extra fields that were not useful.

---

## **WaspLabeler db Tracker 1.1.0**

**(Application: 1.1.0 Templates: 1.1.0 Database: 1.0.8)**

1.1.0 - April 11, 2008

1. Sped up printing and print preview considerably
2. Added a check box to the print drug screen that skips the preview when printing.
3. Updated the 6019 .LAB files to work with the new blackmark media

Scanner Check Digit Note: Many scanners ship with the option to send the check digits for a UPC barcode turned on. If you are using a scanner to input NDC numbers, make sure to use the scanner manual to turn the Send Check Digits option off. When scanning an NDC barcode, if you get extra numbers that do not match the NDC number after you scan you most likely have the check digits turned on.

Template 6019 Blackmark Labels need special setup options to print correctly. Check the website knowledge base or call support for up to date instructions for printing 6019 or any continuous label with blackmarks.

---